



PLEASE FAX/EMAIL COMPLETED FORM TO:

303.242.8474

lab@RockyMtnLabs.com

LAB SUPPLY ORDER FORM

ALLOW UP TO 48 HOURS (2 BUSINESS DAYS) FOR DELIVERY

CLIENT: _____ ORDERED BY: _____

ADDRESS: _____ PHONE #: _____

CITY: _____ STATE: _____ ZIP: _____ DATE: _____

BLOOD COLLECTION SUPPLIES

- Needle Holders, Disposable (250/pk) _____
- Needles, Multisample 21x1 1/4 (48/bx) _____
- Needles, Multisample 22x1 1/4 (48/bx) _____
- Syringe Needles 21x1 (50/bx) _____
- Tourniquets, Non-
Latex (25/bx) _____
- Wipes (200/bx) _____

TUBES

- 5 mL SST (gold) (100/pk) _____
- 6 mL Red Plain Clot (100/pk) _____
- 4 mL Lavender (K-EDTA) (100/pk) _____
- 2.7 mL Lt Blue (Na-Cit Coag) (100/pk) _____
- 7 mL Royal Blue (K-EDTA) (each) _____
- 6 mL Royal Blue (No additive) (each) _____
- 24 Hour Urine Collection Bottles (each) _____

MICROBIOLOGY/ MOLECULAR DIAGNOSTICS

- Diatherix Swab, eSwab (each) _____
- Stool Culture Vial (Orange Lid) (each) _____
- Ova & Parasite w/ Culture Kit (each) _____
- Sterile Stool Collection Cups (each) _____
- GEN-PROBE (Aptima) Urine Kit (CT,GC) (each) _____

CYTOLOGY/PATHOLOGY

- ThinPrep Pap Kit = 25 Vials w/ Brooms _____
- ThinPrep Pap Kit = 25 Vials w/
Brushes/ Spatulas _____
- Tissue Biopsy Vials (Formalin) (each) _____
- RPMI (for Flow Cytometry transport) (each) _____

MISCELLANEOUS

- Glucose Beverage
 - 50 GRAM Lemon-Lime (each) _____
 - 75 GRAM Fruit Punch (each) _____
 - 100 GRAM Orange (each) _____
- Specimen Transport Bags (100/pk) _____
- Sterile Collection Cups (100/pk) _____
- 24 Hour Urine Collection Bottles (each) _____

FORMS

- Test Requisitions _____
- Supply Order Forms _____
- Prenatal Risk Patient Brochures _____
- Fax Form Add-On Tests _____

Additional Requests:

Per Federal Regulations, laboratories are only allowed to provide physician offices with supplies used to collect specimens that will be sent back to that lab for testing. In compliance with those regulations, Rocky Mountain Labs will provide supplies at no charge for the collection of specimens sent to us for testing.